

# Mouse Spermatogenesis-associated serine-rich protein 2 (SPATS2) ELISA Kit

Catalog No: #EK7151

Package Size: #EK7151-1 48T #EK7151-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Spermatogenesis-associated serine-rich protein 2 (SPATS2) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse (Mus musculus)
Other Names	FLJ13117; Nbla00526; P59SCR; SCR59; SPATA10; putative protein product of Nbla00526 serine-rich spermatocytes and round spermatid protein
Accession No.	Q8K1N4
Uniprot	Q8K1N4
GeneID	72572;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

## Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate SPATS2 in samples. An antibody specific for SPATS2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anySPATS2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for SPATS2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of SPATS2 bound in the initial step. The color development is stopped and the intensity of the color is measured.

Note: This product is for in vitro research use only