

CEP55 (Phospho-Ser425) Conjugated Antibody

Catalog No: #C12585



Package Size: #C12585-AF350 100ul #C12585-AF405 100ul #C12585-AF488 100ul

#C12585-AF555 100ul #C12585-AF594 100ul #C12585-AF647 100ul

#C12585-AF680 100ul #C12585-AF750 100ul #C12585-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CEP55 (Phospho-Ser425) Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	CEP55 (Phospho-Ser425) Antibody detects endogenous levels of CEP55 only when phosphorylated at Ser425
Immunogen Description	A synthesized peptide derived from human CEP55 (Phospho-Ser425)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CEP55, Cancer/testis antigen 111, Centrosomal protein of 55 kDa, C10orf3, CT111, Centrosomal protein 55kDa, URCC6, Up-regulated in colon cancer 6
Accession No.	Swiss-Prot#:Q53EZ4NCBI Gene ID:55165NCBI mRNA#:NCBI Protein#:
Uniprot	Q53EZ4
GeneID	55165;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	54
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

Note: This product is for in vitro research use only