

Glutathione peroxidase 7 antibody

Catalog No: #22190



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Glutathione peroxidase 7 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 186 of Human GPX7
Target Name	Glutathione peroxidase 7
Accession No.	Swiss-Prot:Q96SL4Gene ID:2882
Uniprot	Q96SL4
GeneID	2882;
Concentration	0.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 21kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:500

Images



Western Blot analysis of GPX7 expression in transfected 293T cell line by GPX7 polyclonal antibody.
A: Non-transfected lysate.
B: GPX7 transfected lysate.
12% SDS PAGE
Primary antibody diluted at 1: 500



Immunohistochemical analysis of paraffin-embedded H441 xenograft, using GPX7 antibody at 1: 100 dilution.

Note: This product is for in vitro research use only