

PPFIBP2 antibody

Catalog No: #22259

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PPFIBP2 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 120 and 362 of Human PPFIBP2
Target Name	PPFIBP2
Accession No.	Swiss-Prot:Q8ND30Gene ID:8495
Uniprot	Q8ND30
GeneID	8495;
Concentration	0.8mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 98kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images




Sample (30 ug of whole cell lysate)


A: A431

7.5% SDS PAGE

Primary antibody diluted at 1: 1000



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed A431, using PPFIBP2 antibody at 1: 200 dilution.

Background

PPFIBP2, or liprin-beta-2, belongs to the liprin-beta family. See liprin-beta-1 (MIM 603141) for background on liprins.[supplied by OMIM]

Note: This product is for in vitro research use only