

## mOrange Rabbit Polyclonal Antibody

Catalog No: #38079

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	mOrange Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Specificity	The mOrange Rabbit Polyclonal Antibody detects mOrange & mOrange Tag protein.
Target Name	mOrange
Other Names	mOrange
Concentration	1.0mg/ml
Formulation	IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:2000~1:5000

## Images



Western blot analysis of mOrange recombinant protein, using #38079 diluted at 1:5,000.

## Background

mOrange and mOrange2 are extremely bright orange fluorescent protein monomers which can be used as tags or reporters. Both mOrange fluorescent proteins are mutants derived from mRFP1, a monomeric mutant of DsRed.

Note: This product is for in vitro research use only