

Rabbit IgG Isotype Control, PE Conjugated

Catalog No: #I4001P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Rabbit IgG Isotype Control, PE Conjugated |
| Brief Description | Rabbit IgG Isotype Control, PE Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Isotype | IgG |
| Purification | Purified by Protein A. |
| Applications | FCM |
| Conjugates | PE |
| Concentration | 1mg/ml |
| Storage | Store at 4°C for 12 months |

Background

Isotype controls are negative controls that are designed to measure the amount of non-specific background signal caused by the primary antibody. Isotype Control have no specificity for the target protein, yet retain all of the non-specific characteristics of an antibody.

Note: This product is for in vitro research use only