

PAG 1 Antibody

Catalog No: #43478

Package Size: #43748-1 50ug #43748-2 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	PAG 1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB
Species Reactivity	Bo Hu Rt
Specificity	The antibody detects endogenous level of total PAG 1 polyclonal antibody.
Immunogen Description	Recombinant Bovine Pregnancy-associated glycoprotein 1 protein (54-380AA)
Other Names	Pregnancy-associated glycoprotein 1,PAG 1
Accession No.	Swiss-Prot#:Q29432
Uniprot	Q29432
GeneID	281964;
Calculated MW	43KD
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Storage	Store at -20°C

Application Details

WB dilution: □1:500 - 1:5000

Images



Positive WB detected in Raji whole cell lysate, Hela whole cell lysate, Jurkat whole cell lysate

lanes: Pregnancy-associated glycoprotein 1 antibody at 3ug/ml, Secondary Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 43 KDa. Observed band size: 43 KDa.

Product Description

>95%, Protein G purified

Background

Appears to be proteolytically inactive.

Note: This product is for in vitro research use only