

Recombinant Human Galectin-3,LGALS3 (E. coli)

Catalog No: #AP74816



Package Size: #AP74816-1 10ug #AP74816-2 50ug #AP74816-3 500ug #AP74816-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Recombinant Human Galectin-3,LGALS3 (E. coli)
Host Species	E. coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P17931
Uniprot	P17931
GeneID	3958;
Calculated MW	26kDa
Target Sequence	ADNFSLHDALSGSGNPNPQGWPAGWGNQPAGAGGYPGASYPGAYPGQAPPGAYPGQAPPGAYPGAPGAYPGAPAGVYPGPPSGPGAYPSSGQPSATGAYPATGPYGAPAGPLIVPYNLPLPGGVVPRMLITILGTVKPNA NRIALDFQRGNDVAFHFNPRFNENRRVIVCNTKLDNNWGREERQSVFPFESGKPFKIQVLVEPDHFKVAVN DAHLLQYNHRVKKLNEISKLGISGDIDLTSA SYTMI
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, 2mM DTT,pH 7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only