

Recombinant Human IL-36 Receptor Antagonist Protein,IL-36RN,IL-1F5

Catalog No: #AP74818

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP74818-1 10ug #AP74818-2 50ug #AP74818-3 500ug #AP74818-4 1mg

Description

Product Name	Recombinant Human IL-36 Receptor Antagonist Protein,IL-36RN,IL-1F5
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UBH0
Uniprot	Q9UBH0
GeneID	26525;
Calculated MW	16.9kDa
Target Sequence	MVLSGALCFRMKDSALKVLYLHNNQLLAGGLHAGKVIKGEEISVVPNRWLDASLSPVILGVQGGSQCLSCGVG QEPTLTLEPVNIMELYLGAKESKSFTFYRRDMGLTSSFESAAYPGWFLCTVPEADQPVRILTQLPENGGWNAPI TDFYFQQCD
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only