

## Recombinant Human Neurotrophin-3,NT3

Catalog No: #AP74825



Package Size: #AP74825-1 10ug #AP74825-2 50ug #AP74825-3 500ug #AP74825-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Neurotrophin-3,NT3
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P20783
Uniprot	P20783
GeneID	4908;
Calculated MW	13.6kDa
Target Sequence	YAEHKSHRGEYSVCDSESLWVTDKSSAIDIRGHQVTVLGEIKTGNSPVKQYFYETRCKEARPVKNGCRGIDDK HWNSQCKTSQTYVRALTSENNKLVGWRWIRIDTSCVCALSRKIGRT
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 250mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only