

Recombinant Human Sonic Hedgehog,SHH C24II

Catalog No: #AP74830



Package Size: #AP74830-1 10ug #AP74830-2 50ug #AP74830-3 500ug #AP74830-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

Product Name	Recombinant Human Sonic Hedgehog,SHH C24II
Host Species	E.coli
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q15465
Uniprot	Q15465
GeneID	6469;
Calculated MW	19.8kDa
Target Sequence	MIIGPGRGFGKRRHPKKLTPLAYKQFIPNVAEKTLGASGRYEGKISRNSERFKELTPNYNPDIIIFKDEENTGADR LMTQRCKDKLNALAISVMNQWPGVKLRVTEGWDEDGHHSEESLHYEGRAVDITSDRDRSKYGMLARLAVE AGFDWVYYESKAHIHCSVKAENSVAAKSGG
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB, 100mM NaCl, 1mM DTT, pH 7.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only