

Acetyl Lysine Monoclonal Antibody

Catalog No: #HW098



Package Size: #HW098-1 50ul #HW098-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Acetyl Lysine Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC
Specificity	The antibody detects endogenous Acetyl Lysine protein.
Immunogen Type	Protein
Immunogen Description	Purified Protein
Accession No.	Swiss-Prot: NCBI Gene ID:
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1:1000-2000
IHC: 1:200-500

Images



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Acetyl Lysine Monoclonal Antibody.



Western blot analysis of 1) Mouse Brain Tissue, 2) Hela, 3) Hela+TSA Treated using Acetyl Lysine Monoclonal Antibody.

Note: This product is for in vitro research use only