

Gas6 Polyclonal Antibody

Catalog No: #40950



Package Size: #40950-1 50ul #40950-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Gas6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu
Specificity	Gas6 Polyclonal Antibody detects endogenous levels of Gas6 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Gas6.
Target Name	Gas6
Other Names	GAS6; AXLLG; Growth arrest-specific protein 6; GAS-6; AXL receptor tyrosine kinase ligand
Accession No.	Swiss-Prot: Q14393NCBI Gene ID: 2621
Uniprot	Q14393
GeneID	2621;
SDS-PAGE MW	75kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.
Immunohistochemistry: 1/100 - 1/300.
ELISA: 1/40000.
Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using Gas6 Polyclonal Antibody

Note: This product is for in vitro research use only