

HA Tag HRP Antibody

Catalog No: #48052

Package Size: #48052-1 50ul #48052-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HA Tag HRP Antibody
Purification	Immunogen affinity purified
Applications	WB
Immunogen Description	peptide
Conjugates	HRP-conjugated
Other Names	HA epitope tag antibody HA1 antibody HA2 antibody hemagglutinin antibody Hemagglutinin HA1 chain antibody Hemagglutinin HA2 chain antibody
Accession No.	Swiss-Prot#:
Formulation	1*TBS (pH7.4), 0.5%BSA, 40%Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:10,000

Images



Western blot analysis on YPYDVPDYA Tag recombinant protein (55kDa).

Background

Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA molecule corresponding to amino acids 98-106. This antibody is used to detect proteins that are tagged with HA epitope (YPYDVPDYA) and are expressed in prokaryotic and eukaryotic cells.

References

Note: This product is for in vitro research use only