

CD68 Polyclonal Antibody

Catalog No: #41618

Package Size: #41618-1 50ul #41618-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD68 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IHC-p IF
Species Reactivity	Hu Ms Rt
Specificity	CD68 Polyclonal Antibody detects endogenous levels of CD68 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human CD68.
Target Name	CD68
Other Names	CD68; Macrosialin; Gp110; CD68
Accession No.	Swiss-Prot: P34810NCBI Gene ID: 968
Uniprot	P34810
GeneID	968;
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

IHC-p: 1:100-300

Images



Western Blot analysis of mouse spleen cells using CD68 Polyclonal Antibody

Note: This product is for in vitro research use only