

OXSR1(OSR1) antibody

Catalog No: #22738

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	OXSR1(OSR1) antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 251 of OXSR1(OSR1)
Target Name	OXSR1(OSR1)
Accession No.	Swiss-Prot:O95747Gene ID:9943
Uniprot	O95747
GeneID	9943;
Concentration	0.6mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 58kd

Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images




Sample (30 ug of whole cell lysate)


A: 293T

7.5% SDS PAGE

Primary antibody diluted at 1: 10000



Merged with DNA probe



Immunofluorescence analysis of paraformaldehyde-fixed A431, using OXSR1(OSR1) antibody at 1: 200 dilution.

Background

The product of this gene belongs to the Ser/Thr protein kinase family of proteins. It regulates downstream kinases in response to environmental stress, and may play a role in regulating the actin cytoskeleton. [provided by RefSeq]

Note: This product is for in vitro research use only