

Human C1q And Tumor Necrosis Factor Related Protein 3 (C1QTNF3) ELISA Kit

Catalog No: #EK12935

Package Size: #EK12935-1 48T #EK12935-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human C1q And Tumor Necrosis Factor Related Protein 3 (C1QTNF3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Homo sapiens human
Other Names	CORS26; CTRP3; Corcs; Cors; Cartonectin; Collagenous repeat-containing sequence 26 kDa protein
Accession No.	UniProtKB:Q9BXJ4Gene ID:114899
Uniprot	Q9BXJ4
GeneID	114899;
Storage	12 months

Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.127ng/mL

Sample Type:serum, plasma, tissue homogenates or other biological fluids.

Sample Volume: 100 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only