

## Mouse serum

Catalog No: #SR101

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|              |             |
|--------------|-------------|
| Product Name | Mouse serum |
| Storage      | -20°C       |

## Product Description

Source: Normal, healthy, unimmunized mouse

Appearance: Liquid

Filtration: 0.2um

Preservative: None

Stability: Stable for 5 years

Precautions: 1. This product can be kept at 2-8°C for maximum 30 days. Avoid repeated freezing and thawing. If multiple usages are needed, it is recommended to be stored separately in smaller volumes.

2. After taking the serum out of the freezer, it is first placed in a refrigerator at 2-8°C to thaw it, and then at room temperature until it is fully melted. During the melting process, the serum should be gently shaken frequently to melt evenly, however severe vibration should be avoided. 3. The number of serum freeze-thaw cycles should be reduced during use. The serum should not be placed at room temperature or 4 °C refrigerator for prolonged time, as it will affect the serum quality. 4. If there is a flocculent precipitate after the serum is thawed, the serum can be centrifuged at a low speed and the supernatant can be used.

Exogenous materials: No anticoagulant, coagulant or antibiotics were added in the production process.

Total Protein: 55 mg/ml

IgG concentration: 9 mg/ml

pH : 7.80.3

Hemoglobin: 0.4 mg/ml

Note: This product is for in vitro research use only