

CX3CL1(human):Fc(human)(rec.)

Catalog No: #AP79603

Package Size: #AP79603-1 100 ug #AP79603-2 25 ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CX3CL1(human):Fc(human)(rec.)
Brief Description	1545096460
Host Species	Human cell
Purification	>98%
Species Reactivity	Hu
Other Names	Fractalkine; Chemokine (C-X3-C Motif) Ligand 1; Neurotactin; SCYD1
Calculated MW	35 KDa (monomer)
Target Sequence	The extracellular domain of human CX3CL1 (AAH01163.1) (Gln25-Gly100) is fused to the N-terminus of the Fc region of human IgG1.

Product Description

Endotoxin Content:

Note: This product is for in vitro research use only