

Chicken IL-6 ELISA Kit

Catalog No: #EK18717



Package Size: #EK18717-1 48T #EK18717-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chicken IL-6 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Chicken

Application Details

Detect Range:1ng/L -40 ng/L

Sample Type: Serum, plasma, Urine,tissue homogenates and other biological fluids.

Sample Volume:50ul

Assay Time: 3-5h

Detection wavelength:450nm

Note: This product is for in vitro research use only