

AXL (Phospho-Tyr702) Antibody

Catalog No: #12791

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	AXL (Phospho-Tyr702) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	AXL (Phospho-Tyr702) Antibody detects endogenous levels of AXL only when phosphorylated at Tyr702
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human AXL (Phospho-Tyr702)
Target Name	AXL
Modification	Phospho
Other Names	AXL, AXL transforming sequence/gene, Oncogene AXL, AXL oncogene, AXL receptor tyrosine kinase, JTK11, Tyro7, UFO
Accession No.	Swiss-Prot#: P30530NCBI Gene ID: 558
Uniprot	P30530
GeneID	558;
Target Species	human
Calculated MW	98kd
Concentration	1.0mg/ml
Formulation	Rabbit IgG in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Note: This product is for in vitro research use only