

TAS2R41 Antibody

Catalog No: #44656

Package Size: #44656-1 50ul #44656-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAS2R41 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TAS2R41 Antibody detects endogenous levels of total TAS2R41
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TAS2R41
Target Name	TAS2R41
Other Names	TAS2R41, T2R59, T2R41, T2RP27
Accession No.	Swiss-Prot#: P59536NCBI Gene ID:259287
Uniprot	P59536
GeneID	259287;
Calculated MW	36 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only