

FAM9B antibody

Catalog No: #23175



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	FAM9B antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 169 of Human FAM9B
Target Name	FAM9B
Accession No.	Swiss-Prot:Q8IZU0Gene ID:171483
Uniprot	Q8IZU0
GeneID	171483;
Concentration	0.3mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 22kd
Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)
A: HeLa
B: Hep G2
C: Molt-4
D: Raji
12% SDS PAGE
Primary antibody diluted at 1: 5000

Note: This product is for in vitro research use only