

Acetyl-alpha-Tubulin (Lys40) Antibody

Catalog No: #HW169

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Acetyl-alpha-Tubulin (Lys40) Antibody
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Acetyl-alpha-Tubulin (Lys40) Antibody detects endogenous levels of total Acetyl-alpha-Tubulin (Lys40)
Immunogen Type	Peptide
Immunogen Description	Synthesized peptide derived from human alpha-Tubulin around the Acetylation site of Lys40.
Target Name	Acetyl-alpha-Tubulin (Lys40)
Modification	Acetyl
Other Names	TUBA1B, Alpha-tubulin ubiquitous, Alpha tubulin, K-ALPHA-1, Tubulin, alpha 1b, Tubulin alpha-ubiquitous chain, Tubulin, alpha, ubiquitous, Tubulin K-alpha-1, Tubulin alpha chain, Tubulin alpha-1B chain
Accession No.	Swiss-Prot#: P68363NCBI Gene ID:10376
Uniprot	P68363
GeneID	10376;
SDS-PAGE MW	50KD
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis Acetyl-alpha-Tubulin (Lys40) using HeLa whole cell lysates

Note: This product is for in vitro research use only