

p53 (Acetyl-K120) Polyclonal Antibody

Catalog No: #HW186



Package Size: #HW186-1 50ul #HW186-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	p53 (Acetyl-K120) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	This antibody detects endogenous levels of p53 (Acetyl-K120), It doesn't reacte with total protein.
Immunogen Type	protein
Immunogen Description	Synthesized Acetyl peptide derived from human p53. at AA range: K120
Target Name	p53 (Acetyl-K120)
Other Names	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
Accession No.	Swiss-Prot:P04637NCBI Gene ID:7157
Uniprot	P04637
GeneID	7157;
Calculated MW	53 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1B£B'1500-1B£B'12000

Images



Western Blot analysis of 1,mouse-liver 2,hela 3,mouse-brain cells using primary antibody diluted at 1:1000(°C overnight). Secondary antibodyB£B'1Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25B'1£C, 1 hour)

---

Note: This product is for in vitro research use only