

JAK2/3 (Phospho-Tyr966/939) Antibody

Catalog No: #12400



Package Size: #12400-1 50ul #12400-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	JAK2/3 (Phospho-Tyr966/939) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	This antibody detects endogenous levels of JAK2/3 (Phospho-Tyr966/939), It doesn't reacte with total protein.
Immunogen Description	Synthesized Phospho peptide derived from human JAK2/3.at AA range: T966
Target Name	JAK2/3 (Phospho-Tyr966/939)
Other Names	Tyrosine-protein kinase JAK2/JAK3 (EC 2.7.10.2) (Janus kinase 2/Janus kinase 3) (JAK-2/JAK-3)
Accession No.	Swiss-Prot:O60674/P52333NCBI Gene ID:1132/1124
Uniprot	O60674
GeneID	3717;
Calculated MW	125 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1£º500-1£º2000

Images



Western Blot analysis of 1,mouse-liver 2,hela 3,mouse-brain cells using primary antibody diluted at 1:1000(4Λ C overnight).

Note: This product is for in vitro research use only