

## Noggin Polyclonal Antibody

Catalog No: #46870

Package Size: #46870-1 50ul #46870-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Noggin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous Noggin
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 21-70
Target Name	Noggin
Other Names	Noggin
Accession No.	Swiss-Prot:Q13253NCBI Gene ID:9241
Uniprot	Q13253
GeneID	9241;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: IHC 1:50-300

## Images



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

Note: This product is for in vitro research use only