

## COL14A1 Polyclonal Antibody

Catalog No: #46876

Package Size: #46876-1 50ul #46876-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	COL14A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous COL14A1
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 1571-1620
Target Name	COL14A1
Other Names	Collagen alpha-1(XIV) chain (Undulin)
Accession No.	Swiss-Prot:Q05707NCBI Gene ID:7373
Uniprot	Q05707
GeneID	7373;
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Immunohistochemistry: IHC 1:50-300

## Images



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:200

Note: This product is for in vitro research use only