

Recombinant Mouse Zinc finger protein 346(Znf346)

Catalog No: #AP74632




Package Size: #AP74632-1 20ug #AP74632-2 100ug #AP74632-3 1mg

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Zinc finger protein 346(Znf346)
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% as determined by SDS-PAGE.
Immunogen Description	Expression Region:1-294aaSequence Info:Full Length
Other Names	Just another zinc finger protein
Accession No.	Q9R0B7
Uniprot	Q9R0B7
GeneID	26919;
Calculated MW	48.7 kDa
Tag Info	N-terminal 6xHis-SUMO-tagged
Target Sequence	MECPAPDATDAADPGEAGPYKGSEEPEGREPDGVRFDREARRRLWEAVSGAQPAGREEVEHMIQKNQCLF TSTQCKVCCAMLISESQKLAHYQSKKHANKVKRYLAIHGMETIKGDVKRLDSDQKSSRSKDKNHCCPICNMTF SSPAVAQSHYLGKTHAKSLKQKQSTKGAALQQNREMLDPDKFCSLCHSTFNDPAMAQQHYMGKRHRKQE TKLKLMAHYGRLADPAVSDLPAGKGYPCCKIVLNSIEQYQAHVSGFKHKNQSPKTLVTLGSQTPVQTQPTP KDSSTVQD
Formulation	Tris-based buffer50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C,-80°C. The shelf life of lyophilized form is 12 months at -20°C,-80°C.Notes:Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

Binds with low affinity to dsDNA and ssRNA, and with high affinity to dsRNA, with no detectable sequence specificity.

References

"JAZ requires the double-stranded RNA-binding zinc finger motifs for nuclear localization."Yang M., May W.S., Ito T.J. Biol. Chem. 274:27399-27406(1999) Research Topic:Others

Note: This product is for in vitro research use only