

# Mouse anti-Human CD4/CD8/CD3, FITC/PE/PE-Cy5 Conjugated mAb

Catalog No: #28249

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	Mouse anti-Human CD4/CD8/CD3, FITC/PE/PE-Cy5 Conjugated mAb
Host Species	Mouse
Clonality	Monoclonal
Clone No.	4AT4/COV8/4AT3
Isotype	IgG2b/IgG1/IgG2a
Applications	FC
Species Reactivity	Hu
Conjugates	FITC/PE/PE-Cy5
Target Name	CD4/CD8/CD3
Formulation	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL

## Application Details

Vol.per.Test: 20 &mu;l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically add 10ul of this reagent to 100&mu;l of experimental sample with 1 X 10<sup>6</sup> cells per test. Please refer to the detailed protocol when you perform a test.

## Images



Human peripheral blood lymphocytes analyzed with Mouse Anti-Human CD4/CD8/CD3, FITC/PE/PE-Cy5

Note: This product is for in vitro research use only