

## FLJ23560 antibody

Catalog No: #22176

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FLJ23560 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 182 and 413 of Human BBS10
Target Name	FLJ23560
Accession No.	Swiss-Prot:Q8TAM1Gene ID:79738
Uniprot	Q8TAM1
GeneID	79738;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


## Application Details

Predicted MW: 81kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:50-1:500

## Images



Sample (30 ug of whole cell lysate)  
A: Hep G2  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded NCI-N87 xenograft, using BBS10 antibody at 1: 500 dilution.

---

Note: This product is for in vitro research use only