

TMEM8B Antibody

Catalog No: #37186

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TMEM8B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu Ms
Specificity	The antibody detects endogenous levels of total TMEM8B protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human Transmembrane protein 8B
Target Name	TMEM8B
Other Names	NGX6; NAG-5; C9orf127; RP11-112J3.10
Accession No.	Swiss-Prot#: A6NDV4 NCBI Gene ID: 51754Gene Accssion: NP_057530
Uniprot	A6NDV4
GeneID	51754;
Concentration	0.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human colon cancer tissue using #37186 at dilution 1/15.



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using #37186 at dilution 1/15.

Background

May function as a regulator of the EGFR pathway. Probable tumor suppressor which may function in cell growth, proliferation and adhesion.

Note: This product is for in vitro research use only