

POLR3A antibody

Catalog No: #22202

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	POLR3A antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 182 and 452 of POLR3A
Target Name	POLR3A
Accession No.	Swiss-Prot:O14802Gene ID:11128
Uniprot	O14802
GeneID	11128;
Concentration	1.5mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 156kd


Western blotting: 1:500-1:3000

Immunofluorescence: 1:100-1:200

Images



Sample (30 ug of whole cell lysate)
A: HCT116
5% SDS PAGE
Primary antibody diluted at 1: 500



Confocal immunofluorescence analysis (Olympus FV10i) of paraformaldehyde-fixed HeLa, using POLR3A antibody (Green) at 1: 500 dilution and alpha-tubulin antibody (Red) at 1: 2000.

Note: This product is for in vitro research use only