

PEX6 antibody

Catalog No: #22206

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	PEX6 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 180 and 432 of PEX6
Target Name	PEX6
Accession No.	Swiss-Prot:Q13608Gene ID:5190
Uniprot	Q13608
GeneID	5190;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 104kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: A549

7.5% SDS PAGE

Primary antibody diluted at 1: 1000

Background

This gene encodes a member of the AAA (ATPases associated with diverse cellular activities) family of ATPases. This member is a predominantly cytoplasmic protein, which plays a direct role in peroxisomal protein import and is required for PTS1 (peroxisomal targeting signal 1, a C-terminal

tripeptide of the sequence ser-lys-leu) receptor activity. Mutations in this gene cause peroxisome biogenesis disorders of complementation group 4 and complementation group 6. [provided by RefSeq]

Note: This product is for in vitro research use only