

## PLEKHG4 antibody

Catalog No: #22211

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PLEKHG4 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 35 and 226 of PLEKHG4
Target Name	PLEKHG4
Accession No.	Swiss-Prot:Q58EX7 Gene ID:25894
Uniprot	Q58EX7
GeneID	25894;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 131kd

Western blotting: 1:500-1:3000

## Images



Sample (30 ug of whole cell lysate)  
A: JurKat  
5% SDS PAGE  
Primary antibody diluted at 1: 1000

## Background

The protein encoded by this gene contains multiple domains suggestive of a role in intracellular signaling and cytoskeleton dynamics at the Golgi apparatus. Mutations in this gene are associated with spinocerebellar ataxia 16q22-linked. Several alternatively spliced transcript variants, differing

only in the 5' UTR, or encoding a different isoform, have been found for this gene. [provided by RefSeq]

---

Note: This product is for in vitro research use only