

## RAB3C antibody

Catalog No: #22215

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RAB3C antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 197 of RAB3C
Target Name	RAB3C
Accession No.	Swiss-Prot:Q96E17Gene ID:115827
Uniprot	Q96E17
GeneID	115827;
Concentration	0.9mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


## Application Details

Predicted MW: 26kd


Western blotting: 1:500-1:3000

Immunohistochemistry: 1:100-1:500

## Images



Sample (30 ug of whole cell lysate)  
A: H1299  
12% SDS PAGE  
Primary antibody diluted at 1: 1000



Immunohistochemical analysis of paraffin-embedded DLD1 xenograft, using RAB3C antibody at 1: 500 dilution.

## Background

RAB proteins, such as RAB3C, are small GTPases implicated in the regulation of vesicular trafficking between membrane compartments (Cheng et al., 2002 [PubMed 12296628]).[supplied by OMIM]

Note: This product is for in vitro research use only