

Dbs antibody

Catalog No: #22234

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Dbs antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 54 and 307 of Human MCF2L
Target Name	Dbs
Accession No.	Swiss-Prot:O15068Gene ID:23263
Uniprot	O15068
GeneID	23263;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

Predicted MW: 131kd

Western blotting: 1:500-1:3000

Images



Sample (30 ug of whole cell lysate)

A: Molt-4

B: Raji

5% SDS PAGE

Primary antibody diluted at 1: 1000

Note: This product is for in vitro research use only