

## GNPAT Antibody

Catalog No: #37601

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	GNPAT Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total GNPAT protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to residues near the C terminal of human glyceronephosphate O-acyltransferase
Target Name	GNPAT
Other Names	DAPAT; DAP-AT; DHAPAT
Accession No.	Swiss-Prot#: O15228NCBI Gene ID: 8443Gene Accssion: NP_055051
Uniprot	O15228
GeneID	8443;
Concentration	2.5mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol.
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:50-1:200

## Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #37601 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human cervical cancer tissue using #37601 at dilution 1/40.

## Background

This gene encodes an enzyme located in the peroxisomal membrane which is essential to the synthesis of ether phospholipids. Mutations in this gene are associated with rhizomelic chondrodysplasia punctata.

Note: This product is for in vitro research use only