

TC21 antibody

Catalog No: #22241



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	TC21 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 186 of Human RRAS2
Target Name	TC21
Accession No.	Swiss-Prot:P62070Gene ID:22800
Uniprot	P62070
GeneID	22800;
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 23kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:50-1:500

Images



Sample (30 ug of whole cell lysate)
A: A431
12% SDS PAGE
Primary antibody diluted at 1: 5000



Immunohistochemical analysis of paraffin-embedded DLD1 xenograft, using TC21 antibody at 1: 500 dilution.

Note: This product is for in vitro research use only