

PGM3 antibody

Catalog No: #22244



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	PGM3 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 44 and 291 of PGM3
Target Name	PGM3
Accession No.	Swiss-Prot:O95394Gene ID:5238
Uniprot	O95394
GeneID	5238;
Concentration	0.7mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 60kd  
Immunofluorescence: 1:100-1:200


Images



Sample (30 ug of whole cell lysate)  
A: JurKat  
7.5% SDS PAGE  
Primary antibody diluted at 1: 1000



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed HeLa, using PGM3 antibody at 1: 200 dilution.

## Background

Interconverts GlcNAc-6-P and GlcNAc-1-P.

Note: This product is for in vitro research use only