

NAA25 Antibody

Catalog No: #37753

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NAA25 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total NAA25 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human N(alpha)-acetyltransferase 25, NatB auxiliary subunit
Target Name	NAA25
Other Names	NAP1; MDM20; C12orf30
Accession No.	Swiss-Prot#: Q14CX7NCBI Gene ID: 80018Gene Accssion: NP_079229
Uniprot	Q14CX7
GeneID	80018;
Concentration	2.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human thyroid cancer tissue using #37753 at dilution 1/40.



Immunohistochemical analysis of paraffin-embedded Human brain tissue using #37753 at dilution 1/40.

Background

This gene encodes the auxiliary subunit of the heteromeric N-terminal acetyltransferase B complex. This complex acetylates methionine residues that are followed by acidic or asparagine residues.

Note: This product is for in vitro research use only