

CLEC16A Antibody

Catalog No: #37766

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CLEC16A Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total CLEC16A protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human C-type lectin domain family 16, member A
Target Name	CLEC16A
Other Names	Gop-1; KIAA0350
Accession No.	Swiss-Prot#: Q2KHT3NCBI Gene ID: 23274Gene Accssion: NP_056041
Uniprot	Q2KHT3
GeneID	23274;
Concentration	2.9mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #37766 at dilution 1/45.

Background

This gene encodes a member of the C-type lectin domain containing family. Single nucleotide polymorphisms in introns of this gene have been associated with diabetes mellitus, multiple sclerosis and rheumatoid arthritis. Multiple transcript variants encoding different isoforms have been found

for this gene.

Note: This product is for in vitro research use only