

## NME5 antibody

Catalog No: #22253

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	NME5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB
Species Reactivity	Hu Ms
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 188 of NME5
Target Name	NME5
Other Names	NM23H5; NM23-H5
Accession No.	Swiss-Prot:P56597 Gene ID:8382
Uniprot	P56597
GeneID	8382;
Concentration	0.2mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 10% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

## Application Details

Predicted MW: 24kd

Western blotting: 1:500-1:3000

## Images



Sample (30 ug of whole cell lysate)

A: HeLa

12% SDS PAGE

NME5 antibody diluted at 1: 1000

## Background

Does not seem to have NDK kinase activity. Confers protection from cell death by Bax and alters the cellular levels of several antioxidant enzymes

including Gpx5. May play a role in spermiogenesis by increasing the ability of late-stage spermatids to eliminate reactive oxygen species.

---

Note: This product is for in vitro research use only