

SAMD4B Antibody

Catalog No: #37890

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SAMD4B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SAMD4B protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human sterile alpha motif domain containing 4B
Target Name	SAMD4B
Other Names	SMGB; hSmaug2
Accession No.	Swiss-Prot#: Q5PRF9NCBI Gene ID: 55095Gene Accssion: NP_060498/Q5PRF9
Uniprot	Q5PRF9
GeneID	55095;
Concentration	1.2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:10-1:50

Images



Immunohistochemical analysis of paraffin-embedded Human liver cancer tissue using #37890 at dilution 1/20.

Background

SMGB; hSmaug2

Note: This product is for in vitro research use only