

SLC43A2 Antibody

Catalog No: #37946

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SLC43A2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total SLC43A2 protein.
Immunogen Type	Peptide
Immunogen Description	Synthetic peptide corresponding to a region derived from internal residues of human solute carrier family 43 (amino acid system L transporter), member 2
Target Name	SLC43A2
Other Names	LAT4
Accession No.	Swiss-Prot#: Q8N370NCBI Gene ID: 124935Gene Accssion: NP_689559/Q8N370
Uniprot	Q8N370
GeneID	124935;
Concentration	2mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human prostate cancer tissue using #37946 at dilution 1/30.

Background

System L amino acid transporters, such as SLC43A2, mediate sodium-independent transport of bulky neutral amino acids across cell membranes.

Note: This product is for in vitro research use only