

Syntaxin 5 antibody

Catalog No: #22277



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Syntaxin 5 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Purified by antigen-affinity chromatography.
Applications	WB IHC IF
Species Reactivity	Hu
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant protein fragment contain a sequence corresponding to a region within amino acids 1 and 223 of Syntaxin 5
Target Name	Syntaxin 5
Accession No.	Swiss-Prot:Q13190Gene ID:6811
Uniprot	Q13190
GeneID	6811;
Concentration	1mg/ml
Formulation	Supplied in 0.1M Tris-buffered saline with 20% Glycerol (pH7.0). 0.01% Thimerosal was added as a preservative.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.


Application Details

Predicted MW: 30kd
Western blotting: 1:500-1:3000
Immunohistochemistry: 1:100-1:500
Immunofluorescence: 1:100-1:200


Images




Sample (30 ug of whole cell lysate)
A: A431
12% SDS PAGE
Primary antibody diluted at 1: 5000



Immunohistochemical analysis of paraffin-embedded MCF7 xenograft, using Syntaxin 5 antibody at 1: 500 dilution.



Costained with Hoechst 33342



Immunofluorescence analysis of methanol-fixed A431, using Syntaxin 5 antibody at 1: 200 dilution.

Background

Mediates endoplasmic reticulum to Golgi transport.

Note: This product is for in vitro research use only