

Histone H3(Tri-methyl-K4) Mouse Monoclonal Antibody (2E11)

Catalog No: #HW049

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Histone H3(Tri-methyl-K4) Mouse Monoclonal Antibody (2E11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	2E11
Purification	Affinity purification using immunogen.
Applications	WB
Species Reactivity	Hu Rt Ms
Specificity	Antibody detects endogenous Histone H3 (Tri-methyl-K4) protein.
Immunogen Description	A synthetic Tri-methylated peptide corresponding to residues surrounding K4 of human histone H3.
Target Name	Histone H3(Tri-methyl-K4)
Modification	Methyl
Accession No.	Swiss-Prot#:P68431
Uniprot	P68431
GenelD	8350;8351;8352;8353;8354;8355;8356;8357;8358;8968;
SDS-PAGE MW	17kd
Concentration	1.0mg/ml
Formulation	Mouse IgG1 in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500~1:1000

Images



Western blot analysis of 1) Hela, 2) Rat Testis tissue, 3) Raw264.7, using #HW049 diluted at 1:1000.

Background

Histone H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

Note: This product is for in vitro research use only