

Podocin Polyclonal Antibody Cy3 Conjugated

Catalog No: #C06500Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Podocin Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 250-300 383 derived from human Podocin
Conjugates	Cy3
Target Name	Podocin
Other Names	PDCN; SRN1; Podocin; NPHS2
Accession No.	Swiss-Prot#:Q9NP85NCBI Gene ID:7827
Uniprot	Q9NP85
GeneID	7827;
Excitation Emission	512,550nm 570,615nm
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Plays a role in the regulation of glomerular permeability, acting probably as a linker between the plasma membrane and the cytoskeleton.

Note: This product is for in vitro research use only