

Neurotensin Polyclonal Antibody Cy3 Conjugated

Catalog No: #C06517Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Neurotensin Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 14-161 170 derived from human Neurotensin(C-QLYENKPRRPYIL)
Conjugates	Cy3
Target Name	Neurotensin
Other Names	NN; NT; NT N; NTS1; NMN-125; Neurotensin neuromedin N; NTS
Accession No.	Swiss-Prot#:P30990NCBI Gene ID:4922
Uniprot	P30990
GeneID	4922;
Excitation Emission	512,550nm 570,615nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Neurotensin may play an endocrine or paracrine role in the regulation of fat metabolism. It causes contraction of smooth muscle.

Note: This product is for in vitro research use only