

## ECFP Rabbit Polyclonal Antibody

Catalog No: #38077

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ECFP Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purification using immunogen.
Applications	WB
Specificity	The ECFP Antibody detects ECFP & CFP tag protein..
Target Name	ECFP
Other Names	Cyan fluorescent protein; Enhanced Cyan fluorescent protein
Concentration	1.0mg/ml
Formulation	IgG in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:2000~1:5000

## Images



Western blot analysis of ECFP recombinant protein, using #38077 diluted at 1:5,000.

## Background

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Note: This product is for in vitro research use only